




# SwiftX™ ParaBact

## Technical Sheet

► **SwiftX™ ParaBact – designed for bacteria and parasite extraction**

The following components are included in the kit:

- **Buffer BPE** (25mL): Maintains integrity of biological cells and enables efficient capturing of a wide range of bacteria and parasites to our paramagnetic particles.
- **Beads A** (0.75mL): Bifunctional paramagnetic particles enabling species-independent binding of cells as well as removal of inhibitors and cell debris by reverse purification.
- **Buffer BPL** (5mL): Formulated for efficient lysis of cells, particularly in conjunction with the application of heat. The buffer also actively removes RNA from the lysate.
- **Buffer AD** (1mL): Added after cell lysis; renders the DNA lysate compatible with downstream applications such as PCR or isothermal amplification technologies, and nanopore sequencing.



► **Powerful cell lysis**

- Improved cell lysis driven by **alkaline conditions** in addition to heat incubation. Applicable to a wide range of targets, such as various bacteria, parasitic protozoa, fungi, and animal and human host cells
- **Superior** to competitor’s direct lysis methods due to bead-based **inhibitor removal** (reverse purification)

Comparison of lysis efficiencies of Buffer BPL (part of SwiftX™ ParaBact) and Buffer DL (part of SwiftX™ DNA):

Species	Kingdom	Order	Buffer BPL	Buffer DL
<i>Escherichia coli</i>	Bacteria	Enterobacterales	++	++
<i>Staphylococcus aureus</i>	Bacteria	Bacillales	+	o
<i>Pseudomonas aeruginosa</i>	Bacteria	Pseudomonadales	++	++
<i>Mycobacterium tuberculosis</i>	Bacteria	Mycobacteriales	++	++
<i>Streptococcus pneumoniae</i>	Bacteria	Lactobacillales	+	++
<i>Bordetella pertussis</i>	Bacteria	Burkholderiales	++	++
<i>Candida albicans</i>	Fungi	Saccharomycetales	+	++
<i>Aspergillus brasiliensis</i>	Fungi	Eurotiales	++	++
<i>Cryptosporidium parvum</i>	Apicomplexa	Eucoccidiorida	++	+



► **Improved cell capture from liquid specimens**

- The capabilities of SwiftX™ ParaBact for capturing of biological cells have been optimized for **capturing of bacteria and parasites** from various liquid samples types.

Comparison of cell capturing efficiencies of Buffer BPE (part of SwiftX™ ParaBact) and Buffer EN (part of SwiftX™ DNA):

Species	Kingdom	Order	Buffer BPE	Buffer EN
<i>Mycoplasma hyorhinis</i>	Bacteria	Mycoplasmatales	++	++
<i>Staphylococcus aureus</i>	Bacteria	Bacillales	++	++
<i>Listeria monocytogenes</i>	Bacteria	Bacillales	++	+
<i>Streptococcus suis</i>	Bacteria	Lactobacillales	++	-
<i>Streptococcus pneumonia</i>	Bacteria	Lactobacillales	+	o
<i>Leptospira interrogans</i>	Bacteria	Leptospirales	++	++
<i>Borrelia burgdorferi &amp; myamotai</i>	Bacteria	Spirochaetales	++	n.d.
<i>Proteus mirabilis</i>	Bacteria	Enterobacterales	++	++
<i>Klebsiella pneumoniae</i>	Bacteria	Enterobacterales	++	++
<i>Escherichia coli</i>	Bacteria	Enterobacterales	++	o
<i>Salmonella typhimurium</i>	Bacteria	Enterobacterales	++	o
<i>Mycobacterium ulcerans &amp; leprae</i>	Bacteria	Mycobacteriales	++	++
<i>Mycobacterium tuberculosis</i>	Bacteria	Mycobacteriales	++	+
<i>Rhodococcus hoagii</i>	Bacteria	Mycobacteriales	++	+
<i>Anaplasma phagocytophilum</i>	Bacteria	Rickettsiales	++	+
<i>Neoehrlichia spp</i>	Bacteria	Rickettsiales	++	+
<i>Rickettsia spp.</i>	Bacteria	Rickettsiales	++	n.d.
<i>Pseudomonas aeruginosa</i>	Bacteria	Pseudomonadales	++	+
<i>Acinetobacter pittii</i>	Bacteria	Pseudomonadales	++	+
<i>Bordetella pertussis &amp; parapertussis</i>	Bacteria	Burkholderiales	++	++
<i>Aspergillus brasiliensis</i>	Fungi	Eurotiales	++	++
<i>Schistosoma haematobium</i> (eggs)	Animalia	Diplostomida	++	++
<i>Cryptosporidium parvum</i> (oocysts)	Apicomplexa	Eucoccidiorida	+	-
<i>Babesia spp.</i>	Apicomplexa	Piroplasmida	++	o
<i>Leishmania spp.</i>	Euglenozoa	Trypanosomatida	n.d.	++
<i>Plasmodium spp.</i>	Apicomplexa	Haemospororida	n.d.	+
<i>Toxoplasma gondii</i>	Apicomplexa	Eucoccidiorida	n.d.	+

► **Validated for a wide range of sample types**

- ✓ Urine
- ✓ Saliva
- ✓ Sputum
- ✓ Whole blood
- ✓ Cerebrospinal fluid
- ✓ Fine-needle aspirates
- ✓ Wastewater
- ✓ Cell cultures
- ✓ Tissue homogenates
- ✓ Liquid-based cytology media
- ✓ Guanidine-free transport media

