

SwiftX™ DNA Product Sheet



▶ SwiftX DNA enables rapid and easy DNA extraction

- Rapid extraction of DNA from humans, animals as well as from bacteria, protozoa and viruses (lysis as short as 5 minutes)
- Efficient inhibitor removal by magnetic bead-based reverse purification technology
- Enzyme-free DNA extraction
- No spills due to centrifuge-free workflow



▶ Adaptive and easy-to-perform workflow:

- Universal cell concentration step for enhanced sensitivity with liquid samples such as blood, urine, saliva, cell cultures and others.
- Direct extraction by one-step protocol for swabs, tissue, dried blood and other samples




▶ Eco-friendly by reducing waste and time

- No hazardous ingredients, no organic waste – free of guanidine, enzymes, alcohols and phenol
- Minimal worksteps and consumables required
- No refrigeration – storage at ambient temperature



SwiftX™ DNA was designed to match workflow and performance of the discontinued QIAGEN SpeedXtract® Nucleic Acid Kit

SwiftX is a trademark of Xpedite Diagnostics GmbH
Contact information: www.xpedite-dx.com | info@xpedite-dx.com



SpeedXtract® is a trademark of QIAGEN GmbH

