

SwiftX™ Hi-Sense Product Sheet



▶ Rapid and easy extraction of DNA

- Extraction of bacterial, viral and eukaryotic DNA from liquid specimens of human and animal origin in 15 minutes
- Efficient inhibitor removal by magnetic bead-based reverse purification technology
- Centrifugation-free cell and virus capture



▶ Highly sensitive extraction from low-titer samples

- Concentration of analytes and removal of specimen matrix ensures exceptional sensitivity of your assay
- Applicable to liquid samples such as saliva, urine, oral fluid, swabbed samples in transport media and others
- Standard sample volume: 1 mL



▶ Eco-friendly by reducing waste and time

- No hazardous ingredients, no organic waste – free of guanidine, enzymes, alcohols and phenol
- Minimal worksteps and consumables required
- No refrigeration – storage at ambient temperature

