

# SwiftX™ Swabs

## Product sheet



### ▶ SwiftX Swabs is a fast RNA purification kit

- Powerful one-step protocol
- 15 min RNA extraction and inhibitor removal
- From any type of swabs, guanidine-free media (VTM, UTM, saline) and saliva
- CE-IVD certified



### ▶ Highly flexible and easy to use: 1 to 200+ samples per hour

- No investment (CAPEX) in automation required
- Maximizes flexibility in laboratory workflows:
  - no minimum number of samples per extraction run
  - testing can start immediately when samples arrive in lab
- Technicians get trained in minutes due to intuitive protocols



### ▶ Eco-friendly by saving waste, time, consumables and money

- No hazardous ingredients such as guanidine or alcohols
- No organic waste (SwiftX Swabs is free of organic solvents such as ethanol, isopropanol or phenol)
- No dedicated and expensive plastic plates or tubes required
- Ultra-low consumption of tips and tubes (even no pipette tip required if working with swab samples)
- SwiftX Swabs is cheaper than most RNA purification kits



### ▶ Validated for SARS-CoV-2 RNA purification from multiple sample types

- Validated by independent laboratories around the world
- Used worldwide for point-of-care PCR testing and in high-throughput laboratories

