

SwiftX™ Sepsis

REF SXFB-20

Intended use

SwiftX™ Sepsis is intended to be used for extraction and purification of DNA of bacterial and fungal cells as well as leukocytes from up to 10 mL of whole blood. For research use only.

The extracted DNA can be applied to a wide range of analytical methods, such as spectrophotometric and fluorometric analyses, real-time quantitative PCR, and isothermal amplification.

How to use

SwiftX™ Sepsis works very differently than conventional blood extraction products. Therefore, we highly recommend to **have a look at chapter “Read before first use” in the IFU** prior to your actual bench work. This section comprises helpful hints and clues to enhance your user experience and to enable a successful start in the lab. **If something goes wrong, the Troubleshooting guide in the Instructions for Use may be helpful.**

The **Instructions for Use** with all information on SwiftX Sepsis can be downloaded by scanning the QR code, by following the below link, or by sending an email to sales@xpedito-dx.com.



Scan QR code to download the full IFU

Full IFU: <https://t1p.de/7lqaj>



Xpedito Diagnostics GmbH
Lilienthalstr. 2a
85399 Hallbergmoos
Germany

Phone: +49-811-998538-10
E-Mail: info@xpedito-dx.com
Web: www.xpedito-dx.com