



SwiftX™ Virus

Product Sheet



▶ Extraction of viral RNA/DNA from serum and plasma

- Extraction of nucleic acids of tropical viruses in **6 minutes**
- Advanced lysate clearing by **reverse purification**
- **Centrifugation-free** workflow that can be performed anywhere




▶ Simple one-step protocol

- **Direct extraction** of viral RNA and DNA blood serum and plasma
- Virus lysis and **nuclease inactivation** in a single step
- No lengthy purification process required



▶ Eco-friendly by reducing waste

- **Free from hazardous chemicals** such as guanidine, alcohol, phenol
- **Only 2 gram of waste** per extraction



Serum, Plasma
+ Proteinase K
+ Beads A



Extracted
RNA/DNA